



Single use soapy sponge
for a proper cleansing and hygiene of the skin

babybaño® soapy sponges are recommended for daily personal care.
Sponges gently remove the skin dirt and impurities with a pleasant massage effect.

Units per package	Description	Presentation	Units per box	Sponge size	Peso da esponja
1	Polyurethane foam Density 15 Kg/m ³	12 handbags with 200 sponges	2,400 sponges	13 x 7,5 cm ± 5% thickness: 1 cm	1.6 – 2.0 g

TRANSPORT INFORMATION

Box measures: 720 x 585 x 840 mm

Weight box: 8 kg ± 10%

Pallet: 5 boxes

Pallet measures: 120 x 80 x 245 cm

BAR CODE:

GTIN 14: 18426696003528

STORAGE AND SHELF-LIFE

- Store in a cool, dry place out of direct sunlight.
- Expiry date 36 months. Once opened consume within 12 months

PHYSICAL PROPERTIES

The sponges are made of foam polyurethane which provides a soft, effective and lightweight material.

GEL: DESCRIPTION, COMPONENTS AND PROPERTIES

Soapy sponges are impregnated with dermatologic soap. pH 5.5 (± 0,3). Soap tested under pediatric control for infants from 0 to 3 years old.



Especially recommended for children from 0 to 3 years old

INGREDIENTS	FUNCTION
Aqua	Solvent
Sodium Laureth Sulfate	Cleansing, emulsifying, foaming, surfactant
Cocamide DEA	Emulsifying, emulsion stabilising, foam boosting, surfactant, viscosity controlling
Disodium Cocoamphodiacetate	Cleansing, foam boosting, hair and skin conditioning, hydrotope, surfactant
Citric Acid	Buffering, chelating, masking
Parfum	Masking
Lactic Acid, Propylene Glycol, Malic Acid, Rosmarinus Officinalis Extract, Camellia Sinensis Leaf Extract, Citrus Aurantium Amara Flower Extract	Preservative
Sodium Hydroxide	Buffering

QUALIDADE MICROBIOLÓGICA

Limit (Aerobic bacteria, fungi and yeast) per cm ² of product	0 UFC / cm ²
--	-------------------------

Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Candida albicans	Absence
---	---------



INFORMATION

The product complies with the requirements of Regulation (EC) number 1223/2009 of the European Parliament and the Council of 30 November 2009 on cosmetic products.

The skin compatibility of the product has been tested under pediatric and dermatological control. No animal testing.

The production plant is certified with ISO 22716:2007 regarding product documentation and Good Manufacturing Practices for cosmetic products.

The company is certified according to the following standards: UNE-EN ISO 9001:2015, UNE-EN ISO 14001:2015 and UNE 166002:2014, in addition to the Innovative SME seal.

PRECAUÇÕES DE UTILIZAÇÃO

- Manter fora do alcance de crianças.
- Não ingerir.
- Evitar contato com os olhos.
- Uso externo.

MODO DE EMPREGO



1. O gel dermatológico é ativo com um pequeno quantidade de água.



2. Esfregar delicadamente ensaboando a pele.



3. Clarificar com água.



4. Descarte a esponja.

RECOMENDAÇÕES PARA SUA APLICAÇÃO

Gel testado sob controle pediátrico recomendado especialmente para uso em crianças de 0 a 3 anos.

